DIGITAL FIT® TENS/EMS 1200

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User manual

Mode d'emploi

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Elektrostimulationsgerät DIGITAL FIT[®] TENS/EMS 1200

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Electro stimulation unit DIGITAL FIT[®] TENS/EMS 1200

Table of Contents -

Chapter 1 - SAFETY INSTRUCTIONS	17
1.1 Safety precautions	17
1.2 Four safety - contraindications	17
1.3 Precautions.....	18
1.4 Special notices	18
Chapter 2 - PURPOSE OF USE OF DIGITAL FIT[®] TENS/EMS 1200	18
Chapter 3 - OPERATION OF DIGITAL FIT[®] TENS/EMS 1200.....	19
3.1 Operation controls on the front of the device	19
3.2 Operation of DIGITAL FIT [®] TENS/EMS 1200.....	19
3.3 Further adjustments of the device	20
3.3.1 Creating own programs	20
3.3.2 Operation of the therapy memory	21
3.4. Background information of the effectiveness of electrostimulation in the incontinence treatment.....	21
Chapter 4 - PLACEMENT OF ELECTRODES AND PROBES.....	21
4.1 Adjacent placement at the pain area	22
4.2 Placement in the area of the spinal nerves.....	22
4.3 Placement of motor, trigger and acupuncture points	22
4.4 Placement of the electrodes on muscle groups	22
.....4.5 Placement of the probes for the incontinence treatment	22
4.6 Simultaneously usable placement alternatives	22
4.7 Information for attaching the electrodes.....	22
Chapter 5 - SPECIFICATIONS, PROGRAMS, SIGNS AND SYMBOLS.....	23
5.1 Technical dates and further characteristics	23
5.2 Programs of the device	23
5.3 Signs and symbols	25
Chapter 6 - MAINTENANCE, DISPOSAL AND GUARANTEE	25
6.1 Cleaning of the device.....	25
6.2 Replacement of batteries	25
6.3 Maintenance of the Gel-Electrodes.....	26
6.4 Immediate maintenance	26
6.5 Disposal and guarantee	26
Chapter 7 - BASIC EQUIPMENT OF THE DEVICE.....	27
Chapter 8 - EXAMPLES FOR FURTHER PLACEMENTS OF ELECTRODES. 27	
..... 8.1 Electrode placement for pain therapy.....	27
8.2 Examples for further placement of electrodes.....	29

Dear DAVITA® customer,

We are pleased to welcome you as a new customer to our product. You have made the right decision. You have purchased a product which can improve your quality of life for many years.

Chapter 1 – SAFETY INSTRUCTIONS

1.1 Safety precautions

Before the first use of the device the user manual should be read carefully. The following precautions have to be considered:

1. it has to be checked visually if there are any signs of damage;
2. all instructions at the unit have to be considered;
3. all boxes with liquid has to be kept from the device; the device is not protected against the penetration of liquid (IP 20);
4. it has to be checked if the accessories for application of the stimulation (electrodes, connection cable etc.) are available.

1.2 For safety - contraindications

Warning:

The electro stimulation is not allowed to use in the following cases:

1. People with pacemakers, heart diseases or cardiac arrhythmias.
2. In front of the neck because of possible laryngospasm (spasmodic contraction of the muscles of the larynx with narrowing of the vocal apparatus or the glottis);
3. beside the neck because of the sensitivity of the carotid glomus;
4. in the abdominal area and in the sacrolumbal during the pregnancy;
5. in the areas of damaged skin, on mucous membranes, by changes in the local sensitivity, in infectious processes, inflammation, phlebitis and thrombophlebitis;
6. patients with implanted electronic devices (e.g. pacemakers) immediately below the stimulated area;
7. near metallic implants and internal fixators (Prosthesis, osteosynthesis material, spirals, screws, plates), especially while using single phase current (diadynamic, faradic, DC for iontophoresis and galvano therapy);
8. in areas with tumors and in the area of the stomach or the lumbal area with lithiasis (Liver or kidney stones);
9. in the area of dermatitis, tactile hypoesthesia, allergies to substances used during the treatment, thermal anesthesia using iontophoresis.

10. It has to be kept out of the reach of children.

1.3 Precautions

In the following cases it is recommended to consult a doctor and to use the electro stimulator under his observation:

1. Use of the neuromuscular stimulation for pain relief by patients with pain syndroms to determine the cause of disease. Especially if there is no pain relief within a few days after beginning the pain therapy;
2. Application of electrodes on the face;
3. near metallic implants and internal fixators (Prosthesis, osteosynthesis material, spirals, screws, plates);
4. Patients with Parkinson's disease because of possible changes in the sensation level;
5. while suspected or identified pregnancy the stimulation is allowed to be done only in peripheral areas such as legs and arms;
6. Applications of electrodes in areas of tendon and muscle injuries;
7. Patients with spastic paralysis.

During the application of electrodes the following points have to be considered:

1. It is not allowed to put on the eyes.
2. It is not allowed to use directly in sensitive areas such as the middle part of the chest, in groin or armpit area.
3. The heart muscle is not allowed to lie between the plus and minus pole, that means near the heart it is not allowed to put an electrode in front at the upper body and at the same time behind at the back.

1.4 Special notices

1. It is not allowed to use the electro stimulator after a big meal. Because of the increased blood flow into another area of the body the digestive process can be extended.
2. If during the treatment symptoms of tachycardia or extrasystoles should occur or if the stimulation is not tolerated, the treatment should be stopped.
3. It is possible that at the end of the application a slight redness appears at the place where the electrodes were used or that some capillaries rupture. In this case it is recommended to use a lighter program or a lower intensity of stimulation.

Chapter 2 – PURPOSE OF USE OF DIGITAL FIT® TENS/EMS 1200

The electro stimulation device DIGITAL FIT® TENS/EMS 1200 is a very multifunctional easy usable unit, that is especially qualified for the following areas.

1. **Pain therapy:** Acute and chronic pain can be treated effectively with the TENS-programs.


2. **Muscle stimulation:** For the training of the muscles several EMS-programs are available. Muscles can be strengthened and the endurance can be improved.
3. **Massage:** To relax the muscles and to improve the circulation the device has particularly suitable massage programs.
4. **Incontinence treatment:** To treat stress, urge and mixed incontinence special programs are developed. The device can be used successfully for bladder and bowel incontinence.

Chapter 3 – OPERATION OF DIGITAL FIT® TENS/EMS 1200

The electro stimulation device DIGITAL FIT® TENS/EMS 1200 is easy and safe to use with four buttons.

3.1 Operating controls on the front of the device

The following diagram shows the operating controls of the DIGITAL FIT® TENS/EMS 1200. The exact operation of the device is described in this chapter.


3.2 Operation of electro stimulation device DIGITAL FIT® TENS/EMS 1200

1. The cover of the back of the unit has to be opened. 4 pieces of 1.5 volt batteries have to be inserted into the device.
2. To use **programs for pain therapy, muscle stimulation or massage** the two connection cables have to be introduced into the connectors on the top of the device. At each connection cable have to be connected 2 electrodes and placed on the particular place of the skin.

Notes for electrodes placement can be found in Chapter 4 Electrode placement.

If an **incontinence treatment** has to be done a connection cable is inserted into the connector on the top of the device on channel A or B. This cable has to be connected with a probe.

3. The ON/OFF button has to be pressed, to turn on the device.
4. With the two group buttons channel A and B can be operated independently and simply.

5. The ↑ or ↓ arrow of the group button must be pressed repeatedly until the numbers on the left side of the call/mem function (Program memory) start to flash. With the + or – button the desired program for the application can be selected.

The detailed list of the available programs of the device can be find in chapter 5 „Specifications and preset programs“.

The **DIGITAL FIT® TENS / EMS 1200** has 12 TENS-programs for pain therapy, 8 programs for muscle stimulation, 4 massage programs and 4 incontinence programs. In addition 4 programs can be programmed and stored by the user.

6. To adjust the **treatment duration** easily the **time-button** must be pressed so long until the desired treatment duration for example 20 minutes appears below the timer symbol on the top of the display.
7. To adjust the **intensity** the ↑ or ↓ arrow of the group buttons control must be pressed repeatedly until the numbers on the left side of the mA start to flash. The intensity can increased with the + button of the group buttons control and decreased with the – button.
The adjustment of the intensity depends on various individual factors such as skin resistance, skin density e.g.. In all cases, where the stimulation is unpleasant the intensity should be reduced or the treatment has to be finished.
If the intensity was increased to quickly and therefore the application is perceived as too strong the – button can be pressed to stop the increasing of the intensity immediately.
8. The program can be finished at any time by pressing the **ON/OFF-button** and the unit is switched off.

3.3 Further adjustments of the electro stimulation device DIGITAL FIT® TENS/EMS 1200

In addition to the above described settings there is the possibility to create **own programs** and to read and to reset the **therapy memory**.

3.3.1 Creating own programs

1. The channels A and B of the unit are locked (Lock).

To unlock the channels (Unlock), the ↑ or ↓ arrow of the group buttons control must be pressed so often until the digits on the left side of the call/mem function (program memory) start to flash. Then the + button of the group buttons control must be

pressed and hold for 30 seconds. Then the lock beside the channel name will disappear.

2. To create an own program a similar program has to be selected under call/mem. The individual values for mode, width, rate, rise time, time on and rest time can be increased with the + button of the group buttons control or decreased with the – button.
3. To store the new created program the + button of the group buttons control under call/mem has to be pressed and hold for 2 seconds.
4. After creating new programs the channel A and B of the unit should be locked. Therefore the + button of the group buttons control under call/mem has to be pressed and hold for approx. 30 seconds.

Important Notice:

Only the programs 29, 30, 31 and 32 are available for creating own programs. **If more programs are created, existing programs will be overwritten from program 1. This has to be avoided!** In case of overwriting specific programs the parameters of the before existing programs have to be inserted and stored manually.

3.3.2 Operation of the therapy memory

With the therapy memory the cumulative duration of the therapy can be called up and the memory can be reset to 0 again. If the therapy memory shows for example 06 23, the total treatment duration is 62 hours and 30 minutes.

1. To call up the cumulative therapy duration the ↑ or ↓ arrow of the group buttons control has to be pressed so often until the digits on the left side of the call/mem function (program memory) start to flash. Then the two + buttons of the group buttons control have to be pressed and hold for 5 seconds. Below on the display the time for the therapy duration is shown.
2. If the two + buttons of the group buttons control are pressed and hold again simultaneously for 5 seconds, the device will return to the normal operation.
3. In order to reset the therapy memory to 0 again, both – buttons of the group buttons control have to be pressed and hold simultaneously for 5 seconds.

3.4 Background information of the effectiveness of electro stimulation in the in continence treatment

It is assumed that the electro stimulation effects on the nerve pudendi and thus a contraction of the pelvic floor muscles is caused. The sphincter is triggered by electrical current pulse without activity of the patient. In this so called automatic pelvic floor exercises the electro stimulation can cause a muscle hypertrophy, an increase in muscle tone and an improvement of contractility of the pelvic floor and cause a normalization of the reflection pattern of the sphincter.

Chapter 4 – PLACEMENT OF ELECTRODES AND PROBES

To achieve an optimal result with the electro stimulation, it is important to place the electrodes correctly. Because each person reacts differently to electro stimulation and individual needs can differ from the preset settings, the user must find himself the optimal location of the electrode placement.

In the following the recommended electrode placements can be find:

4.1 Adjacent placement at the pain area

With this commonly used placement technique the electrodes are placed closed to or around the localized pain area. One electrode is attached to one side of the pain area. The current flows through or around the pain area. This method is particularly suitable when the pain is in the limb or deep in the tissues.

4.2 Placement in the area of the spinal nerves

These are areas of the body which are traversed by a spinal nerve. With this placement an electrode can be placed at the location of the pain and the other one on the nerve root at the spinal cord. It is stimulated directly in the enervated area.

4.3 Placement of motor, trigger and acupuncture points

A motor point is an entry point of a motor nerve in a muscle. Trigger points can be characterized, that there is a considerable pain by touching these points. Acupuncture points are specific points of the traditional Chinese medicine, which can be stimulated to fight diseases.

These points have a high conductivity tissue and are located at various parts of the body. Normally an electrode is placed directly on the point of pain and the second electrode in the further area of the pain.

4.4 Placement of the electrodes on muscle groups

To perform a muscle training the electrodes must be placed directly to the involved muscle groups, e.g. on the quadriceps or biceps.

4.5 Placement of the probes for the incontinence treatment

The vaginal probe has to be inserted into the vagina.

For women and men the rectal- or mini-vaginal probe has to be introduced into the anus. Also usable for women with a constricted vagina.

4.6 Simultaneously usable placement alternatives

Our electro stimulation device has two independent channels. It is possible to combine several placement alternatives of the electrodes at the same time. Thus for e.g. the first channel can be used for the adjacent placement of electrodes and the second channel can be used for the point therapy.

4.7 Information for attaching the electrodes

1. It should only be used on healthy skin.
2. Before the application the skin has to be washed and dried carefully to achieve a good conductivity.
3. The complete surface of the electrode has to be fixed to the skin. The electrodes should not be used any longer, if they don't adhere to the skin.
4. In the case of skin redness after the electro stimulation treatment in the area under the electrodes an application in the same area has to be avoided if the redness is still existing.
5. The device has to be switched off before attaching or removing the electrode from the skin.

Examples of electrode placement can be found in chapter 8.

Chapter 5 – SPECIFICATIONS, PROGRAMS, SIGNS AND SYMBOLES

5.1 Technical dates and further characteristics

Number channels:	2 separate adjustable channels
Power source:	Battery, 4 x 1.5 volt (not included in delivery)
Wave form:	Biphasic rectangular
Impulse width:	50 – 350 μ S
Intensity:	0 – 60 mA
Impulse frequency:	1 – 200 Hz
Timer:	1 – 99 minutes
Types of applications:	Pain therapy, muscle stimulation, massage and incontinence therapy
Existing Programs:	12 TENS-programs, 8 EMS-programs, 4 massage-programs and 4 incontinence programs for each channel
Storable programs:	4 Programs for each channel
With integrated treatment hour counter and locking of the channels with the adjustments.	
Dimensions:	130 x 70 x 30 mm
Weight without batteries:	115 g
Class II:	Type BF
Protection against anesthetic, flammable gases:	Do not belong to AP or APG protection class
Device for the use of:	Temporary use

5.2 Programs of the DIGITAL FIT[®] TENS/EMS 1200

No.:	Program	Mode	Impulse frequency (Hz)	Impulse wide (μ s)	Time on (Sec.)	Rest time (Sec.)	Rise time (Sec.)
TENS-PROGRAMS							
01	High frequency, conventional TENS programs for immediate pain relief via blocking the pain forwarding	N1	80	200			
02		N1	150	75			

03	Low frequency conventional TENS programs for a longer persistent pain relief	N1	2	250			
04	Burst-program against acute pain	B	2	200			
05	Modulated TENS-programs to prevent a habituation regarding the TENS-therapy which can appear by the conventional programs	N3	100	250			
06		N2	100	150			
07	Conventional TENS-program against arthritis	N1	120	75			
08	Modulated TENS-program against shoulder and neck pain	N3	120	175			
09	Modulated TENS-program against back pain	N4	80	200			
10	TENS-programm against knee pain	N1	2	200			
11	Modulated TENS-Program against karpaltunnel syndrome	N3	150	50			
12	Modulated TENS-program against hip and sciatica pain	N3	120	50			
EMS-PROGRAMS							
13	EMS-program muscle atrophy arms	M	25	200	4	12	1
14	EMS-program muscle atrophy legs	M	35	350	4	12	2
15	EMS-program for mobilisation of wrist, elbow and knees after immobilization	M	35	200	5	6	2
16	EMS-program for reducing the spasticity	N1	2	250			
17	EMS-program for strengthen the arms	M	45	200	6	12	2
18	EMS-program for strengthen the legs	M	45	350	6	12	2
19	EMS-program for endurance	M	20	300	10	10	2
20	EMS-program for muscle relaxation and blood circulation	M	30	350	6	8	1
MASSAGE-PROGRAMS							
21	Massage-program against muscular tension	B	1	250			
22	Deep tissue massage program for improving the local circulation	N4	70	175			
23	Massage-program against muscle sore & cramp, providing oxygen into the muscles	N3	100	150			

24	Lympative massage-program for improving the circulation in laymphatic & venous system	M	20	300	6	6	1
INKONTINENCE-PROGRAMS							
25	URO-program against stress incontinence	M	35	250	4	8	2
26	URO-rogram against Stress incontinence	M	50	200	5	10	2
27	URO-program against urge incontinence	M	10	250	5	3	2
28	URO-program against miexed incontinence	M	20	250	5	10	2
OWN PROGRAMS							
29	Optional storable program						
30	Optional storable program						
31	Optional storable program						
32	Optional storable program						

5.3 Signs and symbols

The following signs and symbols appear in this manual, on the device or on it's accessories:

Caution: Accompanying documents have to bew consulted

CE Mark: Conforms to essential requirements of the Medical Device Directive 93/42/EWG.

Application part of the type BF

DISPOSAL: This product is not allowed to dispose as unsorted municipal waste. It is necessary that such waste has to be collected separately for special treatment. Waste Electrical Electronic Equipment (WEEE)

Chapter 6– MAINTENANCE, DISPOSAL AND GUARANTEE

6.1 Cleaning of the device

For cleaning the device we recommend to use a slightly moistended cloth. Fluids should not be used, because the class of protection (IP20) is not appropriate for this.

6.2 Replacement of batteries

The electro stimulation device is delivered without batteries. It has to be inserted into the device 4 pieces of batteries à 1.5 volt.

A battery indicator will light up on the LCD display when the batteries are almost empty.

If the unit will still be used for a longer time, the power of the batteries decreases further. To avoid the loss of program parameters, the unit can not be operated any longer. Now immediately replace the batteries.

With the ON/OFF-button the stimulator has to be turned off. On the back side of the unit the battery cover has to be removed by pushing the convex and it has to be slided downwards, and the batteries have to be taken out. The new batteries have to be inserted correctly with respect to the polarity (note + and – symbols).

6.3 Maintenance of the Gel-Electrodes

1. For the electrodes maintenance some water has to be applied to the adhesive electrodes after each application. It has to be dried a few seconds before attaching the electrode to the silicone pad. The bag with the electrodes has to be closed carefully to avoid desiccation. It has to be kept at room temperature.
2. It has to be kept in mind that the working life of an electrode depends on correct maintenance as well as on the type of current used and its strength. A clean skin is of great importance. Electrodes which do not stick properly to the skin must be reattached.
3. It has to be applied only to healthy skin. In case of allergies or skin irritation, treatment has to be stopped and a doctor has to be consulted.

6.4 Immediate Maintenance

An immediate maintenance to be carried out by DAVITA® or further qualified staff is required in the following cases:

1. Mechanical strain, e.g. the instrument was dropped and, therefore, damaged;
2. The instrument was overheated after leaving it near a source of heat;
3. Suspicion that water has entered the instrument;
4. Case or other parts of the device are damaged, broken or missing;
5. The unit does not seem to work properly.

6.5 Disposal and guarantee

DAVITA® joins to the system „Grüner Punkt“. All parts of the packaging are recyclable materials and environmentally acceptable to dispose. Batteries and electro devices which are signed with the crossed wheelie bin symbol are not allowed to be thrown into the domestic waste. It exists the obligation to dispose such devices separately. The communal public authority which is responsible for the disposal gives information regarding the possibilities of a regulated disposal.

NOTICE SECURITY TECHNICAL CONTROLS (STK) FOR COMMERCIAL USERS:

According to MP BetriebV § 6 (Safety control) is the commercial user obligated to let conduct a safety control every 2 years. The specific distributor has to be contacted.

DAVITA® guarantees for this device a period of 24 months after date of sale. Wear parts are precluded from the guarantee. The guarantee is valid for material, defects and manufacturing faults and includes the provision of spare parts as well as labour. Guarantee expires if the device has been damaged or used improperly or unauthorized intervention was carried out. The guarantee becomes invalid if the unit has not been used in accor-

dance with the instructions. The purchase receipt as proof has to be kept and it has to be submitted together with the claim for guarantee. DAVITA® will only fulfil the guarantee upon receipt of the purchase proof.

Chapter 7 – BASIC EQUIPMENT OF THE DEVICE

DIGITAL FIT® TENS/EMS 1200, Art.Nr.: 30052

- ✚ 1 piece electro stimulation device
- ✚ 1 piece connection cable set for TENS and TENS/EMS 1200, Art.No.: 31013
- ✚ 1 pouch à 4 pieces self adhesive pads 5 x 5 cm with plug connection for TENS and TENS/EMS 1200, Art.No.: 31011
- ✚ 1 piece packing case
- ✚ 1 piece user manual in English language

Spare parts and accessories

- ✚ 1 piece connection cable set for TENS and TENS/EMS 1200, Art.No.: 31013
- ✚ 10 pouch à 4 pieces self adhesive pads 5 x 5 cm with plug connection for TENS and TENS/EMS 1200, Art.No.: 31012
- ✚ 5 pouch à 4 pieces self adhesive pads 5 x 5 cm with plug connection for TENS and TENS/EMS 1200, Art.No.: 31017
- ✚ 3 pouch à 4 pieces self adhesive pads 5 x 5 cm with plug connection for TENS and TENS/EMS 1200, Art.No.: 31018
- ✚ 1 piece Vaginal applicator, Art.Nr.: 31014 (For the application of incontinence programs necessary) or optional
- ✚ 1 piece Rectal / smaller vaginal applicator, Art.No.: 31015 (For the application of incontinence programs necessary)

Chapter 8– EXAMPLES OF ELECTRODE PLACEMENTS


On the following pages drawings of the electrode placements can be find. These are examples. Further placements are possible.

If 2 electrodes are placed which are controlled via one channel of the device in the particular drawing is spoken of 1. If 4 electrodes are placed which are controlled via two channels of the device in the particular drawing is spoken of 2.

8.1 Electrode placement for pain therapy

Two electrodes have to be connected to the connection cable with the help of plug connections. One male connector of the connection cable has a red colour. This is the (+) pole. The other male connector of the connection cable has a white colour. This is the (-) pole. To achieve optimal results, the electrode which was connected to the white male connector, should be placed near the most painful area.


Arthritis

Application with 2 electrodes (1)

Application with 4 electrodes (2)
Respectively adverse


Shoulder and neck pain

Application with 2 electrodes (1)

Application with 4 electrodes (2)

Back pain (lumbago)


Application with 4 electrodes with horizontal placement

Application with 4 electrodes with vertical placement

Application with 2 electrodes with horizontal placement


Application with 2 electrodes with vertical placement

Knee pains


Application with 4 electrodes

Carpal tunnel syndrome


Application with 4 electrodes

Hip- and sciatic pains

8.2 Examples for further placement of electrodes

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Technical changes can be applied at any time without prior notice. Rev. No. : 11121111, Art. No. : 00470

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DAVITA® Medizinische Produkte GmbH & Co. KG • Postfach 2004 • D-47518 Kleve