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## User manual DAVITON®

Dear DAVITA-customer,

With the DAVITON® you can further improve the effectiveness of the light therapy units.

### Initial operation

On the side of the device there is a **black switch for switching on the unit. To start the DAVITON®, you have to press the red button**, then the control-LED on the front side illuminates and a control tone appears. **The black switch must also be used to terminate the application.**

For a successful light therapy it is necessary that you look directly into the light box every 1–2 minutes for a few seconds. So that you will be reminded the DAVITON® gives after switching on an acoustical signal every 90 seconds for 3 seconds.

Please look directly into the light box for the duration of the acoustical signal. Also between the acoustical signals the eyes should always be opened and there is the possibility, for example to read, to watch television or to do other activities.

### Change of battery

The DAVITON® is working with a commercial 9-volt E-Block battery (not included in delivery). The new processor controlled acoustical system gives up a warning signal in frequent intervals of approx. 0,15 seconds if a status of flat battery is reached. For changing the battery you have to open the battery case on the back side and substitute the old battery and close again the battery case. The DAVITON® is ready to use again. Please regard the legal regulations of the current guideline of disposal for batteries and units.

### Notice

The guarantee is two years from buying date.

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Produkte GmbH & Co. KG