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Dear customers,

Medical explanation of the effect of light therapy

How to use the light therapy units

Some basic information before using the light therapy units

How to carry out light therapy with the light therapy unit LD 72 and LD 100

How to carry out light therapy with the unit PhysioLight LD 220

we are pleased to welcome you as a new customer to our product. You have made the right decision. You have purchased a product which will improve your quality of life for many years.

Today we know that sun light acts like a natural clock for many of our physiological processes. Bright light travels in the form of electric impulses via the eyes and the optical nerves to the nervous system's inner zeitgeber in our brain. This "time giver" processes the information about light and darkness and instructs other parts of the brain to rhythmically distribute chemical messengers. Bright light thus stimulates the production of chemical messengers which influence the mood and suppresses the production of the body's sleep inducing hormone.

The light therapy units LD 72 / LD 110 / LD 220 / LD 440 / LD 880 and LD 1100 can be used for the following symptoms:

1. winter depression
2. non-seasonal depressions
3. premenstrual syndrome („PMS“)
4. sleep disorders

Before use, kindly check that the instrument is in functional order. If there is any sign of damage, immediately contact your reseller. Make sure the ventilation slits are clear of any objects obstructing them. Phototherapy Instruments must only be operated in dry and closed rooms that are free of any explosion risks. The instruments must be free standing and must not be exposed to sun.

The Phototherapy Instruments LD 440, LD 880 and LD 1100 are fixed to a flat wall with the screws and plugs included (LD 440 - distance of 505 mm, LD 880 - distance of 1105 mm and LD 1100 - distance of 1405 mm). Make sure the wall has the appropriate load-bearing capacity. If required, special plugs must be used (in case of cavity walls). Please note that any installation work should be carried out by an electrician.

An open space of at least 20 cm must be left underneath the instrument. Please connect the enclosed power cable into the instrument's plug connector located on the lower right side of the housing. Connect the other end of the cable to an earthed wall socket.

Place the instrument on an even surface and connect the plug with a wall socket. Sit down in front of the instrument keeping a distance of approx. 55 cm / LD 72 and 45 cm / LD 110. The application time is approx. 2 hours / LD 72 and 1 hour / LD 110. Move the tilting switch to „ON“ and the instrument glows. After the application time, switch off the instrument.

Place the instrument on an even surface and put the plug into a wall socket. Sit down in front of the instrument and keep a distance of approx. 50 cm in case of the PhysioLight LD 220. The treatment time with LD 220 is approx. 30 minutes. By turning the control knob of the timer clockwise the Phototherapy instrument begins to glow. After the preset application time has expired, the instrument switches off automatically.

Keep your eyes open during the treatment. You are in a position to read, watch TV or carry out any other activity.

After every 1 - 2 minutes, you should look for a few seconds into the light source. A constant staring into the light source must be avoided. The daily application time can be split into several treatments per day. If you want to cover the Phototherapy instrument after use, allow the instrument to cool down, at least 15 minutes.

Important: In order to be protected against UV-rays, the Phototherapy instrument must only be operated with its original protective filter in place.

How to carry out light therapy with the units PhysioLight LD 440 / LD 880 and LD 1100

One or two person/s can be simultaneously treated with the LD 440 instrument; up to three persons can be simultaneously treated with the LD 880 and up to four persons with the LD 1100 instrument. There must be a distance of approx. 65 cm to the instruments. The application time at this distance and with maximum luminous power is about 30 minutes.

By turning the control knob of the timer clockwise on the instrument's side, the Phototherapy instrument begins to glow. First the selector switch on the instrument's front must be turned into one of the following positions:

Position 1: every 2nd lamp is on

Position 2: all lamps are on

Position 0: no lamp is on, even when the timer is on.

Keep your eyes open during application. You can read or carry out other activities during the treatment. Every 1 or 2 minutes you should look directly into the phototherapy instrument for a few seconds. Constant staring into the instrument should be avoided. After the preset application time has run out, the instrument switches off automatically. It is also possible to switch off the instrument by turning the timer to zero.

Important: the Phototherapy Instruments PhysioLight LD 440, LD 880 and LD 1100 must be operated with the original protective filter in place.

Operation of the light therapy units with Tripod or with the Daviton®

The light therapy unit PhysioLight LD 220 can be used together with a tripod (not included). Also the light therapy units LD 440, LD 880 and LD 1100 can be used in connection with a tripod (not included). The supply of the tripod includes a description on how to instal the PhysioLight LD 220, LD 440, LD 880 and LD 1100 to the tripod. The activated Daviton® (not included) reminds you with an acoustic signal every 90 sec for 3 sec to look into the light source.

Frequency and time of treatment

If the treatment takes place at your practitioner it should initially be carried out 5 times per week for a period of 2 - 3 weeks; thereafter, 3 times per week for a period of 2 - 3 weeks. After that, treatment should take place as needed or when the symptoms occur again. At home, the treatment can take place daily, especially during autumn and winter time. The best treatment time for winter depression is the early morning or the time of the lowest mood. The premenstrual syndrome PMS is treated most effectively in the evening. As for abnormal shifts in a person's day-and-night rhythm (circadian rhythm), i.e. going to bed very early and very early awakening) there are two ways of treating them: if you go to bed very early and get up very early, treatment is most effective in the evening. If you go to bed very late and wake up very late, treatment is best in the morning.

Side eEffects and contraindications, i.e. circumstances that do not allow the use of light therapy units

By using phototherapy the following side effects have been observed: Headaches, nausea, irritability and burning of the eyes can occur for a short period. Caution is recommended for all kinds of eye diseases and an eye specialist should be consulted first. Patients who are being treated with lithium or medicine that can increase the eye's sensitivity towards light, should additionally be supervised by an eye specialist when undergoing phototherapy treatment. An absolute contraindication is not known so far.

The instruments do not need any special maintenance. To clean the outer casing parts of the instrument, we recommend a mild cleansing agent (as used for glass). First, please disconnect the power supply. If lamps must be replaced, this should be carried out by an electrician.

Cleaning, maintenance and recycling

Packaging: DAVITA® partakes in the German recycling system „Grüner Punkt“ (“Green Point”). All packaging parts are made of recycable materials and must be recycled in accordance with the guidelines concerning each material (eco-friendly). The casing of the LD 72 and LD 110 instruments is made of polystyrene (PS), and the front cover is made of polymethylmethacrylate (PMMA). The casings of the Physiolight LD 220 - LD 1100 instruments are made of brushed stainless steel.

Fluorescent lamps as well as instruments marked with the symbol of a crossed out dustbin may not be disposed via the regular garbage. The electrical and electronic components must be recycled separately. Please get in touch with your local authority to find the right way of disposal.

Guarantee

DAVITA® guarantees for this instrument for a period of 24 months after date of sale. The guarantee is valid for material defects and manufacturing faults and includes the provision of spare parts as well as labour. Guarantee expires if the instrument has been damaged or used improperly or unauthorized intervention was carried out. Please keep the purchase receipt as proof and submit it together with you claim for guarantee. DAVITA® will only fulfil the guarantee upon receipt of your purchase proof. Lamps are generally excluded from the guarantee.

Instrument Type:	Light Box LD 72	Light Box LD 110	PhysioLight LD 220	PhysioLight LD 440	PhysioLight LD 880	PhysioLight LD 1100
Article number:	10072	10110	10220	10440	10880	11100
Medical product regarding 93/42/EWG	Yes	Yes	Yes	Yes	Yes	Yes
Simultaneous treatment of:	1 person	1 person	1 person	up to 2 persons	up to 3 persons	up to 4 persons
Housing:	High quality plastic housing					
Housing measures:	H: 61, W: 35, D: 11,5 cm	H: 61, W: 35, D: 11,5 cm	H: 60, W: 40, D: 8 cm	H: 66, W: 70, D: 10 cm	H: 66, W: 130, D: 10 cm	H: 66, W: 160, D: 10 cm
Weight:	3,4 kg	3,9 kg	7,5 kg	20 kg	36 kg	41 kg
Lamps with a light spectrum similar to the sunshine:	2 x 36 watts	2 x 55 watts	4 x 55 watts	8 x 55 watts	16 x 55 watts	20 x 55 watts
Lighting power, distance and therapy length:	10.000 Lux at 20 cm distance for 30 minutes	10.000 Lux at 30 cm distance for 30 minutes	10.000 Lux at 50 cm distance for 30 minutes	10.000 Lux at 65 (40*) cm distance for 30 minutes	10.000 Lux at 65 (40*) cm distance for 30 minutes	10.000 Lux at 65 (40*) cm distance for 30 minutes
Important: Please notice that the light therapy can only successfully used in a distance of at least 45 - 55 cm!	5.000 at 35 cm distance for 60 minutes	5.000 Lux at 45 cm distance for 60 minutes	5.000 Lux at 75 cm distance for 60 minutes	5.000 Lux at 100 (70*) cm distance for 60 minutes	5.000 Lux at 100 (70*) cm distance for 60 minutes	5.000 Lux at 100 (70*) cm distance for 60 minutes
Average lifetime of the lamps regarding the specification of the producer:	2.500 Lux at 55 cm distance for 2 hours	2.500 Lux at 70 cm distance for 2 hours	2.500 Lux at 100 cm distance for 2 hours	2.500 Lux at 150 (100*) cm distance for 2 hours	2.500 Lux at 150 (100*) cm distance for 2 hours	2.500 Lux at 150 (100*) cm distance for 2 hours
Energy costs: (at an average price von 0,12 € / kWh):	approx. 10.000 hours	approx. 10.000 hours	approx. 10.000 hours	approx. 10.000 hours	approx. 10.000 hours	approx. 10.000 hours
Degree of protection:	< 0,02 € per hour	< 0,02 € per hour	approx. 0,02 € per hour	approx. 0,04 € per hour	approx. 0,08 € per hour	approx. 0,10 € per hour
Degree of protection:	IP 20	IP 20	IP 20	IP 20	IP 20	IP 20
Class of protection:	I	I	I	I	I	I
Degree of protection:	B \uparrow Typ B	B \uparrow Typ B	B \uparrow Typ B	B \uparrow Typ B	B \uparrow Typ B	B \uparrow Typ B
Medical device:	class II a	class II a	class II a	class II a	class II a	class II a
Fuse:	1 Fuse 5 x 20 /T 1A L 250V	1 Fuse 5 x 20 /T 2A H 250V	2 Fuse 5 x 20 /T 2A H 250V	2 Fuse 5 x 20 /T 3,15A H 250V	2 Fuse 5 x 20 /T 5A H 250V	2 Fuse 5 x 20 /T 6,3A H 250V
Timer:	Not available	Not available	The units have a timer for stepless presetting of the therapy duration between 0 and 120 minutes.			
Electric ballast units:	Instant starting, flickerfree light, no stroboscopic effect. Buzz free operation. Automatic switch off at the end of lamp live time.					
Light emission:	The selection of lighting source and light filter reduces the energy of light emitted in the Infrared and UV range to a minimum level, so it is harmless to human skin and eyes.					
Storage:	- 10 degrees C up to 50 degrees C and a relative humidity of max. 90 %					
Mains requirements:	230 - 240 V~, 50 - 60 Hz					

*Regarding the LD 440 - LD 1100 the values in brackets are reached by position 1 of the turning knob.

In this position only every second lamp illuminates.

230 - 240 V~, 50 - 60 Hz

approx. 0,08 € per hour

approx. 0,04 € per hour

approx. 10.000 hours

approx. 10.000 hours

approx. 10.000 hours

Technical changes can be applied at any time without prior notice.

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